

Probing Spin Relaxometry in Beesworks[®] Beeswax with Pulsed NMR


Ethan Suresh and Max Misterka

8.13 Experimental Physics I, Spring 2026

Nuclear Magnetism


Electrons, protons, neutrons contain an intrinsic property known as **spin**; superposition of spin up and down

Quantum Representation



Magnetic Moment (J/T): $\mu = \frac{\gamma \hbar}{2}$

Classical Representation



Net Magnetization $\vec{M} = N \langle \vec{\mu} \rangle$

NMR: What happens to the spin population in an external magnetic field?

$$H_{spin} = -\vec{\mu} \cdot \vec{B}$$


Larmor Precession:

$$\omega_L = \gamma B$$

Spin-Lattice Relaxation (T1):

Spin population converges to direction of B


$$M_z(t) = M_0(1 - \exp(-t/T_1))$$



Free Induction Decay:

phase decoherence due to magnetic field inhomogeneity

Spin-Spin Relaxation (T2): neighboring spin-spin dipole interactions exchange energy, causing different local precession rates




$$M_{xy}(t) = M_0 \exp(-t/T_2)$$


Relaxometry of Beeswax: T2 increases with temperature due to decreasing chemical anisotropy

- **Beesworks® Beeswax:** composite organic material made of esters, hydrocarbons, fatty acids (long chains of C, O, H), each with their own relaxation times


Solid (Crystalline) Beeswax
($< 40\text{ C}$)



Semicrystalline Beeswax
($40-60\text{ C}$)




Liquid Beeswax (60-70 C)



Hypothesis: as temperature increases, dipole interaction strength decreases relative to thermal energy and thus chemical environments become more uniform; less phase decoherence

NMR: Using a Spin Echo For Measuring T2, not T2*



Apparatus for pulsed NMR sequences

Transmission Circuitry

- When the switch is open, amplifies input signal and transmits $B_1(t)$ to sample coil

Control Circuitry:

- Generate oscillating RF signal near resonance
- Pulse duration (t_{90} , t_{180}) controlled by pulse programmer


Signal Detection

- Amplify response signal detected by coil
- Superimpose response signal with reference signal


Oscilloscope

- Display beat signal with frequency


$$\omega_L - \omega_I$$



NMR Sample Apparatus



Sources:




Tuning Capacitor

Anatomy of NMR sample and Temperature visualization


Beesworks[®] Beeswax sample



Sample and Tube materials




Thermal Camera visualization



Experimental Procedure for measuring T2 in Beeswax

0. Initial calibration:

PW1, PW2, tuning
capacitor calibrated for a
90-180 spin echo




1. Beeswax Heating

Sample heated to past 70 C
using a hot plate




2. Spin Echo measurements

- Sample temperature T_{cam} measured using thermal camera
- As sample cools from 70 to 45 C, spin echo amplitude is recorded using oscilloscope and phone camera



Repeat steps 1 and 2 for five different delay times τ : 2, 6, 10 16, 18 ms

Example Spin-Echo Signal for measuring T2




Apparatus Settings

- **Tau:**
- **PW1 (t90):** 31 μs $\pm$ 2
- **PW2: (180):** 62 μs $\pm$ 5
- **Averaging:** 16 (necessary to filter out noise)
- Repeat time was not held constant but was long compared to tau (> 180 ms)


Spin Echo Amplitude Measurement

- Peak-to-peak spin echo amplitude obtained digitally by measuring in the recording

Spin-Echo Amplitudes do not fit decaying exponential trend



Quantifying uncertainty in T_{cam} with cooling rate



Numerical Derivative Estimate

$$\frac{dT}{dt}(T_i) = \frac{T_{i+1} - T_{i-1}}{t_{i+1} - t_{i-1}}$$

Characteristic Delay Time Estimate

$$\Delta t \approx 1\text{s}$$

Mean Temperature Uncertainty

$$\sigma_T(T) = \Delta t \cdot \frac{1}{N_\tau} \sum_{k=1}^{N_\tau} \left| \frac{dT}{dt} \right|_k (T)$$

Other systematic errors we didn't account for

1. Time scale (τ) limitations:

- Spin echo disappears after increasing tau past 18 ms; likely due to trigger timing and averaging
- Large systematic fit uncertainty



2. Anomalous saturation effect

- Increases fit uncertainty; fit a possibly due to influence of T_1 (spin-lattice relaxation) and repeat time; $T_1 \sim T_2$

3. Non-equilibrium Cooling

- Cooling is not **quasistatic** ($t_{\text{relax}} \ll t_{\text{cool}}$)
- Multiple phases present at the same time
- Lowers T_{cam} and might lower T_2
- Look into insulated sample holder


4. Apparatus Errors

- Pulse width miscalibration: likely does not change relative trend
- Sample location in holder


5. Temperature Calibration

- $T_{\text{cam}} < T_{\text{real}}$; systematic temperature shift
- Measure T_{real} simultaneously with T_{cam}

Uncertain temperature dependence of T2





Higher temperatures likely have a much lower T_2 ; systematic effect from fit



Are outliers a statistical fluke or physically motivated, eg. A phase transition?

Conclusion: Temperature does seem have an effect on T2 relaxometry in Beesworks® beeswax, but more investigation is required to rigorously determine the nature of this effect.

- **Hypothesis:** T2 increases with temperature due to decreasing chemical anisotropy



- **Results:** Large variation in T2, but inconclusive interpretation
- **Statistical Errors:** Measurement noise, different trials for different tau, fit uncertainty
- **Systematic Error:** Temperature cooling (accounted), low range of tau (unaccounted), etc.

Thank you!

- **Lab partner:** Max Misterka
- **Jlab Instructional Staff:** Prof. Janet Conrad, Jixiang Yang, Simon Rothman, Vinh Tran, Atissa Banuazizi
- **JLab Technical Staff:** Dr. Sean Robinson, Dr. Aaron Pilarcik, Cory Romanov
- **Claude Code** was used to generate code for data analysis
- **ChatGPT** was used for beeswax NMR literature review